


Name

Key

**Microscope Labeling**



**Microscope Use:**

- 15. When focusing a specimen, you should always start with the low power <sup>scanning/</sup> objective.
- 16. When using the high power objective, only the fine adjustment knob should be used.
- 17. The type of microscope used in most science classes is the Compound/Light microscope.
- 18. You should carry the microscope by the base and the arm.
- 19. The objectives are attached to what part of the microscope (it can be rotated to click lenses into place?) Rotating Nose Piece
- 20. A microscope has an ocular objective of 10x and a high power objective of 50x, what is the microscope's total magnification? 500x